# Axio Lab A1 操作手册

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- DIN 58959-4: Quality management in medical microbiology – Part 4: Requirements for investigations using light microscopes
- DIN EN 1335-1: Office work chairs – Part 1: Dimensions – Determination of dimensions
- DIN EN 12464-1: Lighting of work places – Part 1: Indoor work places
- DIN EN 12665: Basic terms and criteria for specifying lighting requirements
- DIN EN 13150: Workbenches for laboratories – Safety requirements and test methods
- DIN EN ISO 15189: Medical laboratories – Particular requirements for quality and competence
- DIN EN ISO 9241-11: Ergonomic requirements for office work with visual display terminals – Guidance on usability
- DIN EN ISO 60601-1-6: Medical-electrical equipment – General requirements for basic safety and essential performance

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- DIN 33402-2: Ergonomics – Human body dimensions – Part 2: Values
- DIN 33406: Workplace dimensions in production
- DIN 33408: Body templates for seats
- DIN 33411: Human physical strength
- DIN 68877: Swiveling work chair - Safety requirements, testing
- DIN EN 614-1: Ergonomic design principles – Part 1: Terms and general principles

- DIN EN 894-1: Ergonomic requirements for the design of displays and control actuators – Part 1: General principles for human interaction with displays and control actuators
- DIN EN 894-3: Ergonomic requirements for the design of displays and control actuators – Part 3: Control actuators
- DIN EN 62079: Preparation of instructions – Structuring, content and presentation
- DIN EN ISO 7250: Basic human body measurements for technological design
- DIN EN ISO 14738: Safety of machinery – Anthropometric requirements for the design of workstations at machinery
- ISO 11226: Ergonomics – Evaluation of static working postures
- SEMI S8-0307: Safety guidelines for ergonomics engineering of semiconductor manufacturing equipment

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





**Binocular ergo photo tube 5-30°/23**  
(100:0/0:100), upright image  
425520-9040-000  
recommended to use:  
**Base plate for microscope stand Axlo Lab,A1**  
430037-9100-000

**Binocular ergo tube 20°/23,**  
reversed image variable,  
continuous vertical adjustment  
425511-0000-000  
mandatory:  
**Base plate for microscope stand Axlo Lab,A1**  
430037-9100-000


**Comfortable binocular ergo photo tube 15°/23 (50:50)**  
vertically adjustable and extendable, upright image  
425520-9050-000  
mandatory:  
**Base plate for microscope stand Axlo Lab,A1**  
430037-9100-000


**Binocular ergo photo tube 20°/23 (100:0/0:100),**  
reversed image variable,  
continuous vertical adjustment  
425512-0000-000  
mandatory:  
**Base plate for microscope stand Axlo Lab,A1**  
430037-9100-000


**Binocular ergo photo tube 20°/23 MAT (100:0/0:100),**  
upright image variable,  
continuous vertical adjustment  
425514-0000-000  
mandatory:  
**Base plate for microscope stand Axlo Lab,A1**  
430037-9100-000


**Binocular tube 30°/23,**  
reversed image  
425520-9000-000

**Binocular tube 20°/23,**  
upright image  
425520-9090-000

**Binocular photo tube Pol 20°/23 (100:0/0:100),**  
upright image  
425520-9100-000  
recommended to use:  
**Base plate for microscope stand Axlo Lab,A1**  
430037-9100-000



**Comfortable binocular ergo tube 8-33°/22,**  
50 mm high, reversed image variable,  
continuous angle adjustment 8-33°  
and vertical adjustment 50 mm  
425522-9040-000  
including:  
**Base plate for microscope stand Axlo Lab,A1**


**Dust cover G**  
459306-0000-000  
**Dust cover for Primo Vert/Axlo Lab,A1**  
415510-1901-000

**Mechanical stage 75x30 R for reflected light**  
432035-9170-000


**Microscope stand Axlo Lab,A1 HAL 35, 4xH, conoscopic, rotary stage 360°**  
430037-9040-000

**Microscope stand Axlo Lab,A1 HAL 35, FL-LED, 5x H**  
430037-9020-000

**Microscope stand Axlo Lab,A1 MAT HAL 50, 5x HD, stage A**  
430037-9050-000

**Fine drive knob with scale, changeable**  
430051-9000-000



**Fine drive knob with scale, changeable**  
430051-9000-000




**Filter holder**  
430037-1037-000

**Interference wide-band filter green,**  
d=32x4  
467803-0000-000  
**Band pass interference filter green 546,**  
d=32x3  
467807-0000-000  
**Heat absorbing filter KG 1, d=32x2**  
467830-0000-000

**Neutral-density filter 0.06, d=32x2**  
467848-0000-000  
**Neutral-density filter 0.25; d=32x2**  
467849-0000-000  
**Polarizing filter 32 mm**  
473600-0000-000

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### 3.1.8 Mounting rotary stage Pol

#### 3.1.8.1 Removing rotary stage pole

- Loosen screw cap (Fig. 3-12/6) from the spring housing (about three rotations).
- Press rotary stage Pol (Fig. 3-12/4) to the front against spring-loaded pin (Fig. 3-12/7), lift it off the stage carrier (Fig. 3-12/5) from the back and remove it upwards.
- Re-tighten screw cap (Fig. 3-12/6).

#### 3.1.8.2 Attaching rotary stage Pol

- I Where necessary, loosen screw cap (Fig. 3-12/6) of spring housing with approx. three rotations.

- I Place rotary stage Pol with the groove of the dovetail (stage bottom) on spring-loaded pin (Fig. 3-12/7).
- I Attach the rotary stage with the clamp screw (Fig. 3-12/8) pointing to the front right.
- I Press the rotary stage Pol to the front against the spring-loaded pin and lower it towards the back into the stage carrier (Fig. 3-12/5), then release it.
- I Re-tighten screw cap (Fig. 3-12/6).

☞ The rotary stage must be mounted so that the vernier scale is on the left side and clamp is on the right side.

#### 3.1.8.3 Dismantling the detachable specimen guide and mounting stage clips

- I Loosen the knurled screw (Fig. 3-12/1) on the specimen guide Pol. Withdraw the specimen guide Pol (Fig. 3-12/2) upwards.
- I Insert the stage clips (Fig. 3-12/9) into the holes provided.


### 3.1.8.4 Removing the stage clips and mounting the detachable specimen guide Pol

- I Remove the stage clips (Fig. 3-12/9) from the rotary stage Pol.
- I Insert the specimen guide Pol (Fig. 3-12/2) with the two cylindrical pins on the underside into the holes provided (Fig. 3-12/3) and tighten the clamp screw (Fig. 3-12/1).

### 3.1.8.5 Centering rotary stage Pol

With high-power objectives centering can be exact only for one selected objective.

All stages are factory-precentered, i.e. while rotating the stage the specimen feature set to the center will remain in the center. If the specimen feature moves off the center of the field of view (Fig. 3-13/5), while rotating the stage, the stage should be re-centered by following this procedure:

- I The KÖHLER illumination on the microscope must be adjusted before centering the stage (see Section 4.1.1).
- I Turn the nosepiece to swing the non-centring objective mount.
- I For centering the stage, use a contrasting specimen and an eyepiece with crossline reticle.
- I Loosen the stage clamping screw (Fig. 3-13/1) and screw cap on the stage carrier (Fig. 3-13/3).
- I Rotate the stage to determine the position of maximum offset of the specimen feature (Fig. 3-13/5, origin of arrow) from the center of the eyepiece reticle.
- I Reset the two centering screws on the stage carrier (Fig. 3-13/2) using a SW 1.5 socket wrench (Fig. 3-13/4) to move the specimen detail by half an arrow length in the direction of the crossline center. Check whether the specimen detail moves when the stage is rotated again; repeat the procedure, when required.

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 The SW 1.5 socket wrenches are located in the storage compartment on the rear side of the microscope stand.

- I When centering is finished, re-tighten screw cap (Fig. 3-13/3).

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### 3.1.8.6 Centring objectives of the polarization stand

The nosepiece 4x Pol is equipped with one fixed and three centerable objective positions. Stage centering of the non-centering objective mount is necessary to ensure that a specimen feature located in the center of the field of view does not drift out while rotating the stage. By centering of the remaining objectives, the specimen feature remains in the center of the field of view even after changing the objective.

- The KÖHLER illumination on the microscope must be adjusted before centering the stage (see Section 4.1.1).

- For centering the stage, use a contrasting specimen and an eyepiece with crossline reticle.
- First turn the nosepiece to swing the noncentering objective mount. Center the rotary stage for the non-centering objective mount as described under 3.1.8.5.
- Turn the nosepiece to move a centering objective mount into the light path.
- Rotate the stage to determine the position of maximum offset of the specimen feature (Fig. 3-14/3, origin of arrow) from the center of the eyepiece reticle.
- Reset the two centering screws on the stage carrier (Fig. 3-14/2) using one SW 1.5 Allen screw-driver (Fig. 3-14/1) each to move the specimen detail by half the arrow length in the direction of the crossline center. Check whether the specimen detail moves when the stage is rotated again; repeat the procedure, when required.
- Center the other two objectives in the same manner

 To maintain this centering accuracy, when replacing it is important not to hold the objective itself, but the knurled ring of the nosepiece to rotate the latter.

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### 3.2.1 Mounting the light intensive co-observer unit

The light-intensive co-observer unit is mounted on the Axio Lab.A1 with a main observer and one or two co-observers in accordance with the separate instructions for use for multi-discussion facilities (Publication No. M60-2-0048).

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### 3.2.2 Mounting polarizer D or filter holder

- I Lift the condenser carrier together with its drive knob upwards as far as it will go.
- I If necessary, unscrew locking and holding pin together with the overview fixture from the condenser carrier.
- I Hold the polarizer or filter holder (Fig. 3-26/4) parallel to the underside of the condenser carrier (Fig. 3-26/1) and screw the holding pin (Fig. 3-26/2) of polarizer (Fig. 3-26/4) with the angled adjusting lever (Fig. 3-26/3) into the front left threaded hole below the condenser carrier as far as it will go.
- I Screw the locking pin (Fig. 3-26/6) with adjusting lever (Fig. 3-26/5) as far as it will go into the rear threaded hole of the condenser carrier.

The other components shown in the system overview at this point must be mounted analogously.

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#### 4.1.4 Setting transmitted light polarization

##### 4.1.4.1 Detecting birefringence

###### (1) Application

The transmitted light polarization method is used for specimens which change the state of polarization of light. These specimens, such as crystals, minerals or polymers, are referred to as birefringent. When these birefringent substances are viewed between crossed polarizers (polarizer  $\perp$  analyzer), they appear bright while their surroundings remain dark.

Birefringent substances are identified in that they show four bright and four dark positions when rotated through  $360^\circ$  between crossed polarizers. Depending on birefringence, thickness and orientation of the specimen, interference colors ranging from gray (mostly with biological specimens) to white, yellow, red and blue appear in this process. These interference colors can be of the first or any higher order.

###### (2) Instrumentation

Polarization methods can be used in the transmitted light on Axio Lab.A1 microscopes for transmitted light polarization and conoscopy.

- Tension-free objectives
- Rotary stage Pol
- Polarizer D (rotatable or fixed)
- Analyzer slide D, fixed, or lambda or lambda/4 compensator
- Depolarizer (for screwing into Axio Lab.A1 tubes) to avoid undesirable polarization effects

The depolarizer is already incorporated in the Axio Lb.A1 stand for conoscopy.

A depolarizer (quartz depolarizer) should be installed in all microscopes used for examining mineral/geological specimens.

A depolarizer suppresses undesirable polarization effects which may occur after the analyzer (e.g. on prism surfaces in the tube), or shifts these to higher orders.

### (3) Setting the microscope

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1 (3)).
- Center rotary stage Pol (Fig. 4-5/1) (see Section 3.1.8.5) and objectives (see Section 3.1.8.6).
- Swivel polarizer (Fig. 4-5/3) into the light path and, if it is rotatable, position it at 0°.
- Insert the analyzer slider (Fig. 4-5/2) into the slit for compensators (if tube does not already have an analyzer). The field of view will appear dark due to the crossed polarizers. With screw-on analyzers, care must be taken to ensure that they are aligned with polarizer D (i.e. crossed position).
- Move the specimen into the field of view and turn it with the rotary stage. As a rule, while being rotated between crossed polarizers, birefringent (anisotropic) specimens will now display the variations in color and intensity described above. However, optically anisotropic substances may also remain dark, if an isotropic direction, e.g. of optically uniaxial or biaxial crystals, is oriented parallel to the direction of observation.


Fig. 4-5 Components for transmitted light polarization

#### 4.1.4.2 Determination of gout and pseudogout

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1 (3)).
- Then swivel the polarizer, rigidly fixed to the lambda plate (445226-0000-000) (Fig. 4-5/3) into the beam path.
- Insert the analyzer slider (Fig. 4-5/2) into the slit for compensators.
- The field of view will appear dark due to the crossed polarizers.
- Swivel the rotatable lambda plate into the beam path and set the metal adjusting lever of the lambda plate to 45°. The gamma direction is orthogonal to the position of the level and is indicated by a white line on the top of the lambda plate.

 The 45° position is at the third white marking on the scale. The scale graduation from one marking to the next is 15°. The 45° position is at the third marking, calculated from the 0° marking. As a further reference point for the correct position of the lever, the Greek letter  $\lambda$  has been engraved on the upper side of the lambda plate, likewise at the 45° position.


Fig. 4-6 Gamma direction

- Select crystals which are oriented in the gamma direction (Fig. 4-6).

#### Evaluation:

If the crystal needles parallel to the gamma direction are yellow and those perpendicular to the gamma direction are blue, they are monosodium urate crystals (gout).

If the crystal needles parallel to the gamma direction are blue and those perpendicular to the gamma direction are yellow, they are calcium pyrophosphate crystals (pseudogout).

 Alternatively, a combination of fixed polarizer (427701-0000-000) and analyzer with fixed lambda plate, 45° (453681-0000-000) can also be used. This offers the advantage that the lambda plate is pre-set to 45°, precluding the risk of incorrect setting. The evaluation is performed in the same way as described above.

#### 4.1.4.3 Determining the direction of oscillation $n_\gamma$

##### (1) Application

The determination of the direction of oscillation of  $n_\gamma$  and  $n_{\gamma'}$  (direction of oscillation with the absolutely or relatively highest refractive index) and  $n_\alpha$  and  $n_{\alpha'}$  (direction of oscillation with the absolutely or relatively lowest refractive index) in relation to the morphological directions, e.g. of crystal surfaces, crystal needles or fibers, provides an important criterion for recognition. It is also employed for the diagnosis of biocrystals (e.g. gout, pseudogout).


Fig. 4-7 Determining the direction of oscillation  $n_\gamma$  using a synthetic fiber as an example

##### (2) Instrumentation

- Eyepiece with crossline reticle
- Tension-free objectives
- Rotary stage Pol (Fig. 4-5/1)
- Polarizer D (rotatable or fixed)
- Screw-in fixed analyzer slide D or  $\lambda/\lambda/4$  compensator combined with analyzer (in Axio Lab.A1 tubes)
- Alignment specimen for polarization microscope (453679-0000-000)

### (3) Setting the microscope

- Set the microscope as for the transmitted light brightfield (see Section 4.1.1), taking care to ensure the correct interpupillary distance in the binocular tube (see Section 3.5.5).
- Center rotary stage Pol (Fig. 4-5/1) and objectives (see Sections 3.1.8.5 and 3.1.8.6).
- Swivel polarizer (Fig. 4-5/3) into the light path and, if it is rotatable, position it at  $0^\circ$ .
- Insert the analyzer slider (Fig. 4-5/2) into the slit for compensators (if tube does not already have an analyzer). The field of view will appear dark due to the crossed polarizers. If not, align the analyzer in the tube or the intermediate plate.
- Set the alignment specimen Pol on the microscope stage and turn to the dark position of the alignment specimen.
- Remove the analyzer and align the crossline reticle with the cracks in the specimen.
- Subsequently reinsert the analyzer and remove the alignment specimen. The pass directions of the polarizer and analyzer will now be parallel to the crossline reticle (polarizer EW, analyzer NS).

 An adjustment of the crossline reticle is not necessary when working with the intermediate plate and the binocular photo tube Pol (425520-9100-000).

- Rotate the rotary stage Pol with the specimen, e.g. a synthetic fiber, until the specimen appears as dark as possible. In this position, the fiber extends parallel to one of the two directions of the crossline reticle.

 Do not change the interpupillary distance on the binocular tube, as the the angle of the crossline reticle to the fiber will be changed.

- Now turn the stage on by  $45^\circ$  so that the longitudinal axis of the fiber is oriented NE-SW (Fig. 4-8). The specimen will display the greatest brightness here (diagonal position). In this position the specimen may have any color.
- Insert the compensator  $\lambda$  (only possible if used with screw-in analyzer in tube or intermediate plate).

Like the specimen, the compensator  $\lambda$  is a birefringent object, albeit with a defined path difference of 550 nm and the principal direction of oscillation  $n_\gamma$  definitely oriented in a NE-SW direction.

By moving compensator  $\lambda$  into the light path, the specimen changes its color. The type of color change depends on the orientation of the specimen (NE-SW or NW-SE).

The changes in color are attributable to optical interference. The interference colors (path differences) in both diagonal positions (NE-SW and NW-SE) of the specimen must be compared in this connection.

The path difference results from the superposition (interference) of the direction of oscillation of the specimen with the direction of oscillation of the compensator  $\lambda$ .

The greater path difference occurs, if the direction of oscillation of the specimen with the absolutely or relatively highest refractive index ( $n_\gamma$  or  $n_{\gamma'}$ ) is parallel to the principal direction of oscillation of the compensator  $\lambda$ . The specimen will then appear, for instance, in greenish-blue (Fig. 4-7/2).

The smallest path difference occurs, if the direction of oscillation of the specimen with the absolutely or relatively lowest refractive index ( $n_\alpha$  or  $n_{\alpha'}$ ) is perpendicular to the direction of oscillation of the compensator  $\lambda$ . The specimen will then appear, for instance, in yellow (Fig. 4-7/3).

#### (4) Conclusions

The grayish-white color appearing first in the bright position in the above example (Fig. 4-7/1) corresponds to a path difference of 150 nm according to the Michel-Lévy color chart (Fig. 4-8).

When the compensator  $\lambda$  is brought into the light path, the non-birefringent "surroundings" of the synthetic fiber appear in a dark red color, which corresponds to the path difference of the compensator of 550 nm (1st order interference color for the path difference of 550 nm corresponds to  $1\lambda$ ).

If the direction of oscillation ( $n_y$  or  $n_y'$ ) of the birefringent specimen to be examined is parallel to the principal direction of oscillation ( $n_y$ ) of the compensator  $\lambda$ , i.e. in NE-SW direction, the path difference of the specimen (e.g. grayish-white: 150 nm) and the path difference of the compensator  $\lambda$  (red: 550 nm). This results in a color change of the specimen from grayish white to greenish-blue (resulting path difference = 700 nm).

If the direction of oscillation of the specimen to be examined is perpendicular to the principal direction of oscillation of the compensator  $\lambda$ , i.e. in NW-SE direction, the path difference of the compensator  $\lambda$  (e.g. grayish-white: 150 nm) is subtracted from the path difference of the compensator (red: 550 nm). In this case, the interference color of the specimen visibly changes from grayish-white to orange (resulting path difference = 400 nm).

 Color charts according to Michel-Lévy are available under Cat. No. 42-312.


Fig. 4-8 Schematic diagram of the color charts according to Michel-Lévy

#### 4.1.4.4 Measuring path differences

The measurement compensators are required for exact measurement of path differences. These return, i.e. compensate, the path difference created by the specimen to zero (black of the first order).

Whereas in the above-described methods the addition or subtraction position was of interest, **solely** the subtraction position is of interest in the measurement.

Path differences in the specimen can assume very small values ( $1/50\lambda$  or 10 nm) and very great values (more than  $10\lambda$  or approx. 5500 nm and more) and with that determine the compensator appropriate for the measurement.

The suitable compensator is determined as follows:

- Set the microscope as for the transmitted light brightfield (see Section 4.1.1), taking care to ensure the correct interpupillary distance in the binocular tube (see Section 3.5.5).
- Accurately position the specimen to be examined on the center of the crosshairs.

- Limit to aperture to a value of about 0.2.
- Turn the rotary stage Pol until the specimen is almost obliterated, i.e. **completely dark**, and set the 45° locking position.
- Rotate the stage **once** (by 45°) so that the specimen is in a diagonal position (bright).

The interference intensity or color leads to the following conclusion:

- If more or less strong interference colors appear on the specimen, the path difference ranges approximately between  $1/2 \lambda$  and  $5 \lambda$ .  
The suitable compensator is:  
**Tilting compensator B 0-5  $\lambda$** .
- If the specimen-side color changes from light gray/white to a strong interference color, when a compensator  $\lambda$  (473704-0000-000) is inserted in the compensator slot, the path difference is  $(1/4 \dots 1/2) \lambda$ .

 A prerequisite for occurrence of the color change effect may be the evaluation in two specimen positions rotated at an angle of 90° from one another, plus a centered stage.

The suitable compensator is:

**Tilting compensator B 0-5  $\lambda$**  or the compensation method according to DE SENARMONT up to  $1 \lambda$  using the **Senarmont compensator 546/4 nm**.

 The compensation method according to DE SENARMONT requires the use of the rotatable analyzer.

- After insertion of the compensator  $\lambda$  and rotation of the object by 90°, the interference color remains white; in this case, however, it is a "higher-order white" and thus the path difference is  $> 5 \lambda$ .  
The suitable compensator is:  
**Tilting compensator K 0-30  $\lambda$**  (Accessory 000000-1115-698)
- A dark gray as appearing interference color suggests very small path differences ( $\lambda/10$  or 54.6 nm).  
The suitable compensator is:  
**Brace-Köhler rotary compensator  $\lambda/10$**  (Accessory 000000-1115-703).

- Insert the compensator into the slit as far as it will go.

The accompanying instructions must be observed for measurement preparation and procedure.

#### 4.1.4.5 Circular polarization contrast

##### (1) Application


Unlike standard polarization contrast, circular polarization contrast does not show any dark (extinction) positions that depend on the angle of rotation (azimuth) of the specimen relative to polarizer or analyzer. This means that when the stage is rotated the same image impression remains, as the light/dark positions are omitted. With optical anisotropy, all transparent specimens display the characteristic interference colors.

##### (2) Instrumentation

- Tension-free objectives
- Rotary stage Pol
- Circular polarizer D (no polarizers may be adapted on the condenser) including corresponding  $\lambda/4$  plate.
- Stationery analyzer slide D or screw-in analyzer (in Axio Lab.A1 tubes).

##### (3) Setting the microscope

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1).
- Center rotary stage Pol or objective (if this has not already been done – see Section 3.1.8.5 or 3.1.8.6).
- Initially do **not** use a specimen for the further settings.
- Swivel in the lower part of circular polarizer D (Fig. 4-9/2) into the beam path until it engages and evaluate the extinction (darkening) of the field of view without the specimen at full light intensity. If this is not optimum, align the analyzer in the tube or intermediate plate as necessary.
- Insert the respective slide 6x20 with  $\lambda/4$  plate (Fig. 4-9/1) as far as it will go into the slot for compensators above the nosepiece.
- Then swivel the upper part of circular polarizer D (Fig. 4-9/4) into the beam path.
- Rotate the lever of the  $\lambda/4$  plate of the circular polarizer D (Fig. 4-9/3) until the extinction is maximum (dark-gray field of view) (lever points 45° to the right).

- 1 Slide 6x20 with  $\lambda/4$  plate
- 2 Lower section of circular polarizer
- 3 Lever for rotating  $\lambda/4$  plate
- 4  $\lambda/4$  plate in upper part of the circular polarizer
- 5 Adjustment slits

**Fig. 4-9 Components for circular polarization contrast**

- An (anisotropic) specimen should not be observed until after the above adjustment.
- Reinsert the specimen to be examined.

In their interference color – dependent on material, specimen thickness and orientation – the specimens appear constant and independent of stage rotation.

For a high-contrast image with higher objective enlargements (from approx. 20x) the illumination aperture must be reduced to a value between 0.15 and 0.20, i.e. the aperture diaphragm must be closed accordingly.

The effect of the  $\lambda/4$  plate (Fig. 4-9/4) can be undone by either swiveling it out of the light path or turning it with lever (Fig. 4-9/3) into one of its two click-stop positions.

#### **4.1.5 Setting transmitted light polarization with the conoscopy stand**

#### **4.1.6 Determining the optical character of crystals**

For the classification (and thus identification) of crystalline material – instead of the observation of the specimen itself – the analysis of an interference image in the objective pupil provides the more valuable information. This image is visible in the eyepiece when an additional lens (so-called Bertrand lens) is switched on. Alternatively, the auxiliary microscope or a diopter may be used to view the interference image.

In contrast to orthoscopy, this is referred to as conoscopy, because the illumination is ideally provided by a wide open cone. In practice this means that the aperture diaphragm is fully open and the objective should likewise have a high aperture.

##### **(1) Application**

The crystal diagnosis is for determining the optical character of transparent and weakly absorbent crystals. This method is also referred to as conoscopy.

Its main application is the classical mineral microscopy. However, synthetic crystals, industrial minerals and plastics (e.g. films) can also be identified and characterized.

##### **(2) Instrumentation**

Conoscopic viewing is preferably carried out on the Axio Lab.A1 microscope for transmitted light conoscopy.

- Tension-free objectives; recommended:  
N-Achroplan 50x/0.8 Pol objective or  
EC Plan-Neofluar 40x/0.9 Pol objective
- Rotary stage Pol
- Polarizer D (rotatable or fixed)
- 0.9 Pol condenser

##### **(3) Setting the microscope for conoscopy**

In the case of uniaxial crystals, the most favorable orientation for conoscopic viewing is obtained with those specimen features (e.g. of a thin section) that in orthoscopic viewing change the brightness as little as possible. In this case, the direction of viewing and the optical axis are parallel. The same applies to biaxial crystals if viewed in or approximately in the direction of one of the two optical axes.


Fig. 4-10 Axio Lab.A1 for transmitted light conoscopy

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1).
- Place the specimen on the stage and focus on it.
- Swivel the analyzer into the beam path (**on** position) with rotary knob **A** (Fig. 4-10/2). The direction of oscillation can be changed using the setting wheel (Fig. 4-10/4) of the analyzer.

#### CAUTION

The movements of rotary knobs **A** and **BL** and the respective setting wheels are coupled with one another. Only **one** control element should therefore be operated at a time and the movement of the other should not be inhibited or blocked. Mechanical damage may otherwise occur.

If rotary knob **BL** is set to the **on** position, if it is not already at the **on** position rotary knob **A** is automatically carried.

If, on the other hand, rotary knob **A** is set to the **off** position, if it is not already at the **off** position rotary knob **BL** is automatically carried.

- Place a selected crystal in the center of the crossline reticle.
- Swivel in objective N-Achroplan 50x/0.8 Pol or EC Plan-Neofluar 40x/0.9 Pol and focus with the focusing drive.
- If necessary, close the luminous-field aperture to avoid superimposition of the axial figure by axial figures of neighboring crystals. The smallest crystal range that can be faded out is 170  $\mu\text{m}$ .
- Switch on Bertrand lens **BL** (Fig. 4-10/1) (Position **on**). The axial figure will appear in the field of view.
- Bring the axial figure into focus with setting wheel (Fig. 4-10/5).

#### (4) Evaluation

Crystalline anisotropic specimens can be separated into optical uni- and biaxial, in each case with "optically positive" or "negative" character.

**Uniaxial** crystals display a **black cross** when the optical axis is parallel to the direction of view. Depending on the size of the birefringence and specimen thickness, concentrically arranged colored **interference rings** (so-called isochromes) may appear (see also Fig. 4-11 second row).

This cross remains closed when the stage is rotated. Depending on the section it may lie within or outside the displayed objective pupil.

With **optically biaxial** crystals, the cross resolves into two **dark hyperbola branches** (the so-called isogyres) depending on stage rotation, which are surrounded by colored interference patterns depending on the amount of birefringence and specimen thickness (suggestive of the figure "8").

Inserting a compensator  $\lambda$  (473704-0000-000) or  $\lambda/4$  (473714-0000-000) or a wedge compensator  $0-4 \lambda$  (000000-1140-663) in the compensator slot with the initial state of the axial figure being as illustrated in Fig. 4-11 results in the following changes in color shown schematically (blue and yellow areas) to the axial figure, thus allowing differentiation in "optically positive" and "optically negative".

	Optically uniaxial		Optically biaxial		
	Positive	Negative	Positive	Negative	
$\lambda$ -Plate (white $\rightarrow$ blue $\rightarrow$ yellow)					+ = blue - = yellow
Quartz wedge (Direction of motion at insertion)					Direction of movement 
$\lambda/4$ plate (position of black spots)					

**Fig. 4-11 Determining optical character**

In the case of less favourable sections in which the cross-hair center is optically uniaxial or the isogyres are optically biaxial specimens outside the objective pupil, an assessment is possible as follows:

- If the black isogyres are **straight** and they run parallel to the pupil (in relation to the cross-hairs), the specimen is **optically uniaxial**.
- If the black isogyres are **curved lines** which wander on a circular path through the pupil, the specimen is **optically biaxial**.

Paying appropriate attention, also such axial figures not illustrated in Fig. 4-11 can be interpreted.

Axial figures can often be better displayed with circular polarization. Particularly the axis angle of optically biaxial specimens (quasi distance between the isogyres) can be determined more clearly. The optical character can also be determined. For this purpose the compensator  $\lambda$  (6 x 20), arranged in the compensator slot, is used.

Two storage compartments for 6x20 slides are located on the reverse side of the conoscopy stand.

#### 4.1.6.1 Demonstrating birefringence with the Axio Lab for conoscopy

##### (1) Application

The transmitted light polarization method is used for specimens which change the state of polarization of light. These specimens, such as crystals, minerals or polymers, are referred to as birefringent. When these birefringent substances are viewed between crossed polarizers (polarizer  $\perp$  analyzer), they appear bright while their surroundings remain dark.

Birefringent substances are identified in that they show four bright and four dark positions when rotated through  $360^\circ$  between crossed polarizers. Depending on birefringence, thickness and orientation of the specimen, interference colors ranging from gray (mostly with biological specimens) to white, yellow, red and blue appear in this process. These interference colors can be of the first or any higher order.

##### (2) Instrumentation

On the Axio Lab.A1 microscope for transmitted light conoscopy:

- Tension-free objectives
- Rotary stage Pol
- Polarizer D (rotatable or fixed)
- Compensator  $\lambda$  or  $\lambda/4$

The depolarizer is already incorporated in the Axio Lb.A1 stand for conoscopy.

A depolarizer (quartz depolarizer) should be installed in all microscopes used for examining mineral/geological specimens.

A depolarizer suppresses undesirable polarization effects which may occur after the analyzer (e.g. on prism surfaces in the tube), or shifts these to higher orders.

### (3) Setting the microscope

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1 (3)).
- Center rotary stage Pol (Fig. 4-12/1) (see Section 3.1.8.5) and objectives (see Section 3.1.8.6).
- Swivel polarizer (Fig. 4-12/3) into the light path and, if it is rotatable, position it at  $0^\circ$ .
- Swivel the analyzer into the beam path and adjust on the setting wheel until the field of view is dark.(Fig. 4-12/2)
- Move the specimen into the field of view and turn it with the rotary stage. As a rule, while being rotated between crossed polarizers, birefringent (anisotropic) specimens will now display the variations in color and intensity described above. However, optically anisotropic substances may also remain dark, if an isotropic direction, e.g. of optically uniaxial or biaxial crystals, is oriented parallel to the direction of observation.


Fig. 4-12 Components for transmitted light polarization on conoscopy stand

#### 4.1.6.2 Determination of gout and pseudogout

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1 (3)).
- Swivel the polarizer, rigidly fixed or rotatable (Fig. 4-12/3) into the beam path. Set the rotatable polarizer to the  $0^\circ$  position.
- With the Axio Lab for polarization push the analyzer (453681-0000-000) into the 6x20 slider mount. (Fig. 4-12/2)
- With the AxioLab for conoscopy swivel the analyzer into the beam path and bring it into a crossed position using the setting wheel. Additionally, insert compensator 6x20 (473704-0000-000) into the slider mount 6x20.
- The field of view will appear dark due to the crossed polarizers.
- Select crystals which are oriented in the gamma direction (Fig. 4-13).


Fig. 4-13 Gamma direction

**Evaluation:**


If the crystal needles parallel to the gamma direction are yellow and those perpendicular to the gamma direction are blue, they are monosodium urate crystals (gout).

If the crystal needles parallel to the gamma direction are blue and those perpendicular to the gamma direction are yellow, they are calcium pyrophosphate crystals (pseudogout).

**4.1.6.3 Determining the direction of oscillation  $n_\gamma$**

**(1) Application**

The determination of the direction of oscillation of  $n_\gamma$  and  $n_{\gamma'}$  (direction of oscillation with the absolutely or relatively highest refractive index) and  $n_\alpha$  and  $n_{\alpha'}$  (direction of oscillation with the absolutely or relatively lowest refractive index) in relation to the morphological directions, e.g. of crystal surfaces, crystal needles or fibers, provides an important criterion for recognition. It is also employed for the diagnosis of biocrystals (e.g. gout, pseudogout).

**Fig. 4-14 Determining the direction of oscillation  $n_\gamma$  using a synthetic fiber as an example**

**(2) Equipment configuration for Axio Lab for conoscopy**

- Eyepiece with crossline reticle
- Tension-free objectives
- Rotary stage Pol (Fig. 4-12/1)
- Polarizer D (rotatable or fixed)
- Compensator lambda or lambda/4 as required
- Alignment specimen for polarization microscope (453679-0000-000)

### (3) Setting the microscope

- Set the microscope as for the transmitted light brightfield (see Section 4.1.1), taking care to ensure the correct interpupillary distance in the binocular tube (see Section 3.5.5).
- Center rotary stage Pol (Fig. 4-5/1) and objectives (see Sections 3.1.8.5 and 3.1.8.6).
- Swivel polarizer (Fig. 4-5/3) into the light path and, if it is rotatable, position it at  $0^\circ$ .
- Swivel the analyzer into the beam path and bring it into a crossed position using the setting wheel (Fig. 4-5/2). The field of view will appear dark due to the crossed polarizers.
- Set the alignment specimen Pol on the microscope stage and turn to the dark position of the alignment specimen.
- Swivel out the analyzer and align the crossline reticle with the cracks in the specimen.
- Subsequently swivel the analyzer back in and remove the alignment specimen. The pass directions of the polarizer and analyzer will now be parallel to the crossline reticle (polarizer EW, analyzer NS).

 An adjustment of the crossline reticle is not necessary when working with the intermediate plate and the binocular photo tube Pol (425520-9100-000).

- Rotate the rotary stage Pol with the specimen, e.g. a synthetic fiber, until the specimen appears as dark as possible. In this position, the fiber extends parallel to one of the two directions of the crossline reticle.

 Do not change the interpupillary distance on the binocular tube, as the angle of the crossline reticle to the fiber will be changed.

- Now turn the stage on by  $45^\circ$  so that the longitudinal axis of the fiber is oriented NE-SW (Fig. 4-15). The specimen will display the greatest brightness here (diagonal position). In this position the specimen may have any color.
- Inserting compensator  $\lambda$  (473704-0000-000).

Like the specimen, the compensator  $\lambda$  is a birefringent object, albeit with a defined path difference of 550 nm and the principal direction of oscillation  $n_\gamma$  definitely oriented in a NE-SW direction.

By moving compensator  $\lambda$  into the light path, the specimen changes its color. The type of color change depends on the orientation of the specimen (NE-SW or NW-SE).

The changes in color are attributable to optical interference. The interference colors (path differences) in both diagonal positions (NE-SW and NW-SE) of the specimen must be compared in this connection.

The path difference results from the superposition (interference) of the direction of oscillation of the specimen with the direction of oscillation of the compensator  $\lambda$ .

The greater path difference occurs, if the direction of oscillation of the specimen with the absolutely or relatively highest refractive index ( $n_\gamma$  or  $n_{\gamma'}$ ) is parallel to the principal direction of oscillation of the compensator  $\lambda$ . The specimen will then appear, for instance, in greenish-blue (Fig. 4-14/2).

The smallest path difference occurs, if the direction of oscillation of the specimen with the absolutely or relatively lowest refractive index ( $n_\alpha$  or  $n_{\alpha'}$ ) is perpendicular to the direction of oscillation of the compensator  $\lambda$ . The specimen will then appear, for instance, in yellow (Fig. 4-14/3).


Fig. 4-15 Schematic diagram of the color charts according to Michel-Lévy

#### (4) Conclusions

The grayish-white color appearing first in the bright position in the above example (Fig. 4-14/1) corresponds to a path difference of 150 nm according to the Michel-Lévy color chart (Fig. 4-15).

When the compensator  $\lambda$  is brought into the light path, the non-birefringent "surroundings" of the synthetic fiber appear in a dark red color, which corresponds to the path difference of the compensator of 550 nm (1st order interference color for the path difference of 550 nm corresponds to  $1 \lambda$ ).

If the direction of oscillation ( $n_y$  or  $n_x$ ) of the birefringent specimen to be examined is parallel to the principal direction of oscillation ( $n_y$ ) of the compensator  $\lambda$ , i.e. in NE-SW direction, the path difference of the specimen (e.g. grayish-white: 150 nm) and the path difference of the compensator  $\lambda$  (red: 550 nm). This results in a color change of the specimen from grayish white to greenish-blue (resulting path difference = 700 nm).

If the direction of oscillation of the specimen to be examined is perpendicular to the principal direction of oscillation of the compensator  $\lambda$ , i.e. in NW-SE direction, the path difference of the compensator  $\lambda$  (e.g. grayish-white: 150 nm) is subtracted from the path difference of the compensator (red: 550 nm). In this case, the interference color of the specimen visibly changes from grayish-white to orange (resulting path difference = 400 nm).

 Color charts according to Michel-Lévy are available under Cat. No. 42-312.

#### 4.1.6.4 Measuring path differences with the Axio Lab for conoscopy

The measurement compensators are required for exact measurement. These return, i.e. compensate, the path difference created by the specimen to zero (black of the first order).

Whereas in the above-described methods the addition or subtraction position was of interest, **solely** the subtraction position is of interest in the measurement.

Path differences in the specimen can assume very small values ( $1/50 \lambda$  or 10 nm) and very great values (more than  $10 \lambda$  or approx. 5500 nm and more) and with that determine the compensator appropriate for the measurement.

The suitable compensator is determined as follows:

- Set the microscope as for the transmitted light brightfield (see Section 4.1.1), taking care to ensure the correct interpupillary distance in the binocular tube (see Section 3.5.5).
- Accurately position the specimen to be examined on the center of the crosshairs.

- Limit to aperture to a value of about 0.2.
- Turn the rotary stage Pol until the specimen is almost obliterated, i.e. **completely dark**.
- Rotate the stage **once** (by 45°) so that the specimen is in a diagonal position (bright).

The interference intensity or color leads to the following conclusion:

- If more or less strong interference colors appear on the specimen, the path difference ranges approximately between  $1/2 \lambda$  and  $5 \lambda$ .  
The suitable compensator is:  
**tilting compensator B 0-5  $\lambda$** .
- If the specimen-side color changes from light gray/white to a strong interference color, when a compensator  $\lambda$  (473704-0000-000) is inserted in the compensator slot, the path difference is  $(1/4 \dots 1/2) \lambda$ .

 A prerequisite for occurrence of the color change effect may be the evaluation in two specimen positions rotated at an angle of 90° from one another, plus a centered stage.

The suitable compensator is:

**tilting compensator B 0-5  $\lambda$**  or the compensation method according to DE SENARMONT up to  $1 \lambda$  using the **Senarmont compensator 546/4 nm**.

 The compensation method according to DE SENARMONT requires the use of the rotatable analyzer.

- After insertion of the compensator  $\lambda$  and rotation of the specimen by 90°, the interference color remains white; in this case, however, it is a "higher-order white" and thus the path difference is  $> 5 \lambda$ .

The suitable compensator is:

**tilting compensator K 0-30  $\lambda$**  (Accessory 000000-1115-698)

- A dark gray as appearing interference color suggests very small path differences ( $\lambda/10$  or 54.6 nm).  
The suitable compensator is:  
**Brace-Köhler rotary compensator  $\lambda/10$**  (Accessory 000000-1115-703).

- Insert the compensator into the slit as far as it will go.

The accompanying instructions must be observed for measurement preparation and procedure.

#### 4.1.6.5 Circular polarization contrast with Axio Lab for conoscopy

##### (1) Application

Unlike standard polarization contrast, circular polarization contrast does not show any dark (extinction) positions that depend on the angle of rotation (azimuth) of the specimen relative to polarizer or analyzer. This means that when the stage is rotated the same image impression remains, as the light/dark positions are omitted. With optical anisotropy, all transparent specimens display the characteristic interference colors.


##### (2) Instrumentation

- Tension-free objectives
- Rotary stage Pol
- Circular polarizer D (no polarizers may be adapted on the condenser) including corresponding  $\lambda/4$  plate.

##### (3) Setting the microscope

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1).
- Center rotary stage Pol or objective (if this has not already been done – see Section 3.1.8.5 or 3.1.8.6).
- Initially do **not** use a specimen for the further settings.
- Swivel the analyzer into the beam path.
- Swivel in the lower part of circular polarizer D (Fig. 4-16/2) into the beam path until it engages and evaluate the extinction (darkening) of the field of view without the specimen at full light intensity. If this is not optimum, align the analyzer as necessary.
- Insert the respective slide 6x20 with  $\lambda/4$  plate (Fig. 4-16/1) as far as it will go into the slot for compensators above the nosepiece.
- Then swivel the upper part of circular polarizer D (Fig. 4-16/4) into the beam path.
- Rotate the lever of the  $\lambda/4$  plate of the circular polarizer D (Fig. 4-16/3) until the extinction is maximum (dark-gray field of view) (lever points 45° to the right).

- 1 Slide 6x20 with  $\lambda/4$  plate
- 2 Lower section of circular polarizer
- 3 Lever for rotating  $\lambda/4$  plate
- 4  $\lambda/4$  plate in upper part of the circular polarizer
- 5 Adjustment slits

**Fig. 4-16 Components for circular polarization contrast on conoscopy stand**

- An (anisotropic) specimen should not be observed until after the above adjustment.
- Reinsert the specimen to be examined.

In their interference color – dependent on material, specimen thickness and orientation – the specimens appear constant and independent of stage rotation.

For a high-contrast image with higher objective enlargements (from approx. 20x) the illumination aperture must be reduced to a value between 0.15 and 0.20, i.e. the aperture diaphragm must be closed accordingly.

The effect of the  $\lambda/4$  plate (Fig. 4-16/4) can be undone by either swiveling it out of the light path or turning it with lever (Fig. 4-16/3) into one of its two click-stop positions.

#### **4.1.7 Setting transmitted light polarization for conoscopic observation – determining the optical character of crystals**

For the classification (and thus identification) of crystalline material – instead of the observation of the specimen itself – the analysis of an interference image in the objective pupil provides the more valuable information. This image is visible in the eyepiece when an additional lens (so-called Bertrand lens) is switched on. Alternatively, the auxiliary microscope or a diopter may be used to view the interference image.

In contrast to orthoscopy, this is referred to as conoscopy, because the illumination is ideally provided by a wide open cone. In practice this means that the aperture diaphragm is fully open and the objective should likewise have a high aperture.

##### **4.1.7.1 Application**

The crystal diagnosis is for determining the optical character of transparent and weakly absorbent crystals. This method is also referred to as conoscopy.

Its main application is the classical mineral microscopy. However, synthetic crystals, industrial minerals and plastics (e.g. films) can also be identified and characterized.

##### **(1) Instrumentation**

Conoscopic viewing is preferably carried out on the Axio Lab.A1 microscope for transmitted light conoscopy.

- Tension-free objectives; recommended:  
N-Achroplan 50x/0.8 Pol objective or  
EC Plan-Neofluar 40x/0.9 Pol objective
- Rotary stage Pol
- Polarizer D (rotatable or fixed)
- 0.9 Pol condenser

##### **(2) Setting the microscope for conoscopy**

In the case of uniaxial crystals, the most favorable orientation for conoscopic viewing is obtained with those specimen features (e.g. of a thin section) that in orthoscopic viewing change the brightness as little as possible. In this case, the direction of viewing and the optical axis are parallel. The same applies to biaxial crystals if viewed in or approximately in the direction of one of the two optical axes.

- Set the microscope as in the transmitted light brightfield according to KÖHLER (see Section 4.1.1).
- Swivel polarizer (Fig. 4-12/3) into the light path and, if it is rotatable, position it at 0°.
- Swivel the analyzer into the beam path and bring into a crossed position with the setting wheel. (The field of view will now appear dark)
- Place the specimen on the stage and focus on it.
- Swivel the analyzer into the beam path (**on** position) with rotary knob **A** (Fig. 4-17/2). The direction of oscillation can be changed using the setting wheel (Fig. 4-17/4) of the analyzer.


Fig. 4-17 Axio Lab.A1 for transmitted light conoscopy

#### CAUTION

The movements of rotary knobs **A** and **BL** and the respective setting wheels are coupled with one another. Only **one** control element should therefore be operated at a time and the movement of the other should not be inhibited or blocked. Mechanical damage may otherwise occur.

If rotary knob **BL** is set at the **on** position, rotary knob **A** is automatically carried if it is not already in the **on** position.

If, on the other hand, rotary knob **A** is set to the **off** position, if it is not already at the **off** position rotary knob **BL** is automatically carried.

- Place a selected crystal in the center of the crossline reticle.
- Swivel in objective N-Achroplan 50x/0.8 Pol or EC Plan-Neofluar 40x/0.9 Pol and focus with the focusing drive.
- If necessary, close the luminous-field aperture to avoid superimposition of the axial figure by axial figures of neighboring crystals. The smallest crystal range that can be faded out is approx. 170 µm.
- Switch on Bertrand lens **BL** (Fig. 4-17/1) (Position **on**). The axial figure will appear in the field of view.
- Bring the axial figure into focus with setting wheel (Fig. 4-17/5).

#### 4.1.7.2 Evaluation

Crystalline anisotropic specimens can be separated into optical uni- and biaxial, in each case with "optically positive" or "negative" character.

**Uniaxial** crystals display a **black cross** when the optical axis is parallel to the direction of view. Depending on the size of the birefringence and specimen thickness, concentrically arranged colored **interference rings** (so-called isochromes) may appear (see also Fig. 4-11 second row).

This cross remains closed when the stage is rotated. Depending on the section it may lie within or outside the displayed objective pupil.

With **optically biaxial** crystals, the cross resolves into two **dark hyperbola branches** (the so-called isogyres) **depending on stage rotation**, which are surrounded by colored interference patterns depending on the amount of birefringence and specimen thickness (suggestive of the figure "8").

Inserting a compensator  $\lambda$  (473704-0000-000) or  $\lambda/4$  (473714-0000-000) or a wedge compensator  $0-4 \lambda$  (000000-1140-663) in the compensator slot with the initial state of the axial figure being as illustrated in Fig. 4-18 results in the following changes in color shown schematically (blue and yellow areas) to the axial figure, thus allowing differentiation in "optically positive" and "optically negative".

	Optically uniaxial		Optically biaxial		
	Positive	Negative	Positive	Negative	
$\lambda$ -Plate (white $\rightarrow$ blue $\rightarrow$ yellow)					+ = blue - = yellow
Quartz wedge (Direction of motion at insertion)					Direction of movement 
$\lambda/4$ plate (position of black spots)					

**Fig. 4-18 Determining the optical character**

In the case of less favorable sections in which the cross-hair center is optically uniaxial or the isogyres are optically biaxial specimens outside the objective pupil, an assessment is possible as follows:

- If the black isogyres are **straight** and they run parallel to the pupil (in relation to the cross-hairs), the specimen is **optically uniaxial**.
- If the black isogyres are **curved lines** which wander on a circular path through the pupil, the specimen is **optically biaxial**.

Paying appropriate attention, also such axial figures not illustrated in Fig. 4-18 can be interpreted.

Axial figures can often be better displayed with circular polarization. Particularly the axis angle of optically biaxial specimens (quasi distance between the isogyres) can be determined more clearly. The optical character can also be determined. For this purpose the compensator  $\lambda$  (6 x 20), arranged in the compensator slot, is used.

Two storage compartments for 6x20 slides are located on the reverse side of the conoscopy stand.

## 4.2 Lighting and contrasting method in reflected light

### 4.2.1 Adjusting the reflected light brightfield according to KÖHLER

#### (1) Application

Reflected light brightfield microscopy is the simplest and most common optical microscopy method for examining opaque samples or specimens, e.g. material sections or wafers.

For a true-to-object imaging, indirect ray bundles, i.e. ray bundles diffracted and scattered on the specimen details, are of major importance in addition to the so-called direct ray bundles. The higher the portion of indirect bundles of rays (aperture), the more realistic the microscopic image according to ABBE will be.

The bundled light emitted by the reflected light unit is reflected on a color-neutral beam splitter and then passes through the objective which focuses the beam on the surface of the sample (so-called condenser function). The objective collects the light reflected by the specimen and together with the tube lens it generates the microscopic intermediate image which can then be visually observed or objectively documented.

#### (2) Instrumentation

Reflected light brightfield viewing is only possible with the stand for reflected light.

- Reflector module brightfield ACR P&C for reflected light in the reflector turret

#### (3) Setting the reflected light brightfield


- The microscope has been started up correctly as described in Section 3.
- The microscope is switched on.
- Adjust the light intensity by turning the regulator (Fig. 4-19/4).
- Position a high-contrast reflected light specimen on the microscope stage.
- Swivel in the 10x objective on the nosepiece (Fig. 4-19/3).
- On the nosepiece (Fig. 4-19/6) swivel in the position with the reflector module brightfield.
- Bring the specimen into focus with the focusing drive (Fig. 4-19/5). If possible, always focus away from the specimen in order to avoid a collision between the objective and the specimen.
- Set the knurled wheel of aperture diaphragm **A** (Fig. 4-19/1) in the middle position (about half open/closed).
- Adjust (reduce) the knurled wheel of luminous-field aperture **F** (Fig. 4-19/2) until the luminous-field aperture is visible in the field of vision.
- Use the focusing drive to adjust the focus on the edge of the luminous-field aperture.
- Now open the luminous-field aperture until the latter just disappears behind the edge of the field of view.
- For aperture diaphragm adjustment (image contrast) remove an eyepiece from the tube socket and look into the socket with the naked eye, or use the auxiliary microscope instead of the eyepiece. This functions only with sufficiently reflecting specimens.

- For specimens with medium contrast characteristics, set the aperture diaphragm with the knurled wheel (Fig. 4-19/1) to about 2/3 to 4/5 of the exit pupil diameter of the objective.

In most applications, this aperture diaphragm setting provides optimum contrast at almost ideal resolution, and is therefore the best compromise for the human eye.

- Then reinsert the eyepiece, adjust the focus with the coaxial coarse and fine focusing knob and adjust the brightness to the reflected light specimen.

 Never use the aperture diaphragm for controlling image brightness. Use the adjusting knob (Fig. 4-19/4) for illumination intensity!

- 1 Knurled wheel of aperture diaphragm A
- 2 Knurled wheel of luminous-field diaphragm F
- 3 Nosepiece
- 4 Light intensity control
- 5 Focusing drive
- 6 Reflector turret

**Fig. 4-19** Microscope settings in reflected light brightfield

## 4.2.2 Adjusting the reflected light darkfield

### (1) Application

The reflected light darkfield method is used for examining incompletely reflecting surfaces with different degrees of reflectivity (ideal reflected light brightfield specimens), i.e. with scratches, ruptures, pores or other disruptions to the even surface. All these light-scattering details light up brightly in the darkfield, whereas the even surface remains dark.

### (2) Instrumentation

Observations in the reflected light darkfield can only be made on Axio Lab.A1 microscopes for reflected light.

- Epiplan-Neofluar, EC Epiplan-Neofluar, Epiplan objectives with the additional designation “HD”
- Reflector module darkfield ACR P&C for reflected light

 The stand for reflected light is equipped with a built-in darkfield stop.

### (3) Adjusting the reflected light darkfield

- Adjust the microscope as described in Section 4.2 for the reflected light brightfield. In order to avoid reflexes, the displayed luminous-field aperture should be located slightly beyond the edge of the field of view.
- If used, remove the 6x20 compensator slide.
- Swivel in the objective position with darkfield objective (HD) on the nosepiece.
- If necessary, swivel in the reflector module darkfield on the reflector turret.
- Completely open the aperture diaphragm and switch off or remove the neutral filter as necessary.
- Place the specimen on the stage and sharpen the image.

### 4.2.3 Adjusting reflected light polarization – Proof of bireflexion and reflexion-pleochroism

#### (1) Application

Reflected light polarization is a further contrasting method for cut surfaces of mineral ore, coal, ceramic products, certain metals and alloys, as depending on the orientation of the crystals and object details, the cut surfaces often react differently when reflected in linear polarized light

The illumination light is linear polarized by the polarizer before passing through the objective onto the specimen surface, where it is reflected. Here the beam parts experience phase differences depending on the structure and polarization optical rotations which, when passing through the analyzer, are displayed in different shades of grey. The grey can be converted into a color contrast with the aid of a compensator with  $\lambda$ -plate.

With objectives of very low magnification, a rotatable  $\lambda/4$  plate arranged in front of the objective (Antiflex cap) permits the reflections to be eliminated even with "dark" specimen surfaces, which otherwise would be unavoidable.

#### (2) Instrumentation

Observations in the reflected light darkfield can only be made on Axio Lab.A1 microscopes for reflected light.

- Rotary stage Pol
- Epiplan-Neofluar Pol, EC Epiplan-Neofluar Pol, Epiplan Pol objectives
- C DIC/DIC/TIC ACR P&C or DIC/Pol ACR P&C or DIC Red I ACR P&C reflector module or Pol ACR P&C reflector module in reflector turret
- Analyzer slide D, fixed or compensator Lambda, 6x20 or Lambda/4, 6x20

#### (3) Adjusting reflected light polarization

- Adjust the microscope as described in Section 4.2 for the reflected light brightfield.
- Swivel reflector module P&C (for DIC or Pol) on the reflector turret into the beam path and insert the analyzer slide (or lambda, lambda/4 compensator) into the 6x20 compartment.
- Insert a specimen, adjust the desired magnification level, focus and observe the specimen in the polarization contrast now present while turning the rotary stage Pol.

A specimen is birefractive when the details display differences in brightness and color which change when the stage is rotated.

For samples with low bireflexion it is advisable to use the analyzer with lambda plate, rotatable.

Pleochroism can be detected when the color of the specimen changes when the stage is rotated (reflected light polarizer turned on, analyzer turned off).

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### C

Cable holder.....	26, 28, 30, 32, 34
Care .....	107
Carrying handle.....	26, 28, 34
Centering screw for condenser .....	26, 28, 30, 34
Checking the instrument .....	108
Circular polarization contrast .....	85, 96
Coarse adjustment .....	26, 28, 30, 32, 34
Color charts.....	83, 94
Color glass holder.....	77
Comfortable ergo tube.....	30
Components .....	37
Optional .....	62
Standard.....	46
Condenser.....	26, 28, 30, 34, 43, 57, 64, 72, 75
Condenser carrier .....	57
Conoscopy .....	34, 87, 98
Contrasting techniques.....	14
Control and functional elements.....	26, 28, 30, 32, 34, 36, 37
Co-observer unit.....	62

### D

Darkfield .....	75, 103
Darkfield stop.....	75
Depolarizer.....	78
Description of the instrument .....	14
Determining crystal character.....	87, 98
Dimensions.....	21
Diopter.....	49

Direction of oscillation .....	81, 92
Drive length.....	52
<b>E</b>	
Ergo photo tube .....	38
Ergo tube .....	36, 38, 39, 40
Ergonomics.....	12, 36
Eyecups .....	50
Eyeglass protection ring .....	50
Eyepiece reticle .....	49, 71
Eyepieces.....	26, 28, 30, 32, 34, 49
<b>F</b>	
Filter holder .....	62
Filter mount.....	42
Filter slide .....	45
Filter slide reflected light .....	32
Fine adjustment .....	26, 28, 30, 32, 34
FL/TL toggle switch .....	65
Fluorescence.....	30, 105
Focusing drive.....	26, 28, 30, 32, 34, 74
Friction torque .....	53
Front lens .....	57, 72
<b>G</b>	
Gout.....	80, 91
<b>H</b>	
Halogen lamp 12 V 50 W .....	60
Height stop on condenser carrier .....	74
<b>I</b>	
Industrial property rights.....	118
Instrument care .....	107
Instrument maintenance .....	108
Instrument safety.....	7
Intended use.....	14
Interpupillary distance.....	38, 70
<b>K</b>	
KÖHLER.....	72, 75
<b>L</b>	
LED lamp .....	58
LED module .....	61
Light intensity .....	26, 28, 30, 32
Light intensity control .....	34
Light sources .....	22
Lighting and contrasting method .....	72, 101
Luminous-field diaphragm.....	26, 28, 30, 32, 34, 73, 77
<b>M</b>	

Mains switch .....	65
Maintenance .....	108
Mechanical stage.....	26, 30, 32, 36, 41, 52
Microscope adjustment .....	68
Microscope stages .....	41
Modulator disk .....	64
<b>N</b>	
Nosepiece.....	26, 28, 30, 32, 34, 45, 50, 56, 72
<b>O</b>	
Objective centering.....	56
Objectives.....	50, 55
ON/OFF switch.....	26, 28, 30, 32, 34
Operation.....	72
Overview fixture .....	44, 63
<b>P</b>	
Phase contrast .....	76
Phase stop.....	77
Photo tube .....	26, 28, 32, 34, 37, 48
Polarization .....	28, 78, 87, 98, 104
Polarizer .....	44, 62, 78, 79, 80, 82, 90, 91, 93
Power cord.....	47
Power supply.....	65
Pseudogout .....	80, 91
Push&click module .....	51
<b>R</b>	
Reflected light .....	30, 32, 101, 103, 104, 105
Reflected light brightfield .....	101
Reflected light darkfield .....	103
Reflected light fluorescence .....	105
Reflected light illumination .....	30, 32
Reflected light polarization .....	104
Reflector turret .....	30, 32, 45, 51
Reflexion-pleochroism .....	104
Replacing the fuses.....	108
Reticle .....	49
Rotary stage Pol.....	28, 34, 42, 54, 55, 78, 81, 85, 90, 92, 96
<b>S</b>	
Service.....	112
Setting up .....	46
Slot .....	26, 28, 30, 32, 34
Specimen guide.....	42, 54
Specimen holder.....	41, 42
Specimen holding plate .....	41
Stage carrier .....	26, 28, 30, 32, 34, 52, 54, 55, 57

Stage clip .....	54
Stand .....	26, 28, 30, 32, 34, 36
Stand models .....	26
Stand, ergonomic .....	36
START-UP .....	46
Switching on/off.....	65
System overview.....	16
<b>ANNEX</b>	
Axio Lab.A1 Index Carl Zeiss	
M60-2-0058 e 04/2010 430037-7044-001 117	
<b>T</b>	
Technical data .....	21
Toggle switch FL/TL.....	30
Tool.....	46
Tool flap .....	26, 28, 30, 32, 34
Transmitted light .....	26, 28, 30, 34, 72, 75, 76, 78, 87, 98
Transmitted light brightfield.....	72
Transmitted light illumination .....	30, 34
Transmitted light lamp.....	26, 28
Transmitted light phase contrast .....	76
Transmitted light polarization .....	78, 90
Troubleshooting .....	109
Tube.....	26, 28, 30, 32, 34, 37, 49
TÜV certificate .....	67
<b>U</b>	
Unpacking.....	46
<b>V</b>	
Vertical adjustment.....	30, 40
Vertical adjustment of condenser.....	26, 28, 30, 34
Vertical adjustment of tube.....	39
Viewing height .....	23, 38, 70
<b>W</b>	
Warning label .....	10
Warranty .....	13
Weight .....	21
White light LED lamp .....	58
Workspace configuration.....	66

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© 除非明确授权，否则不允许传播，复制该文件，使用或通告其内容。违反该手册须做出赔偿。
### **Carl Zeiss MicroImaging GmbH**